



DNA sequencing applications require the use of four or more photomultiplier tubes to detect the fluorescence from DNA fragments labeled with fluorescent dyes — each dye indicating the presence of a DNA fragment with one of the four DNA bases (T, A, G, C). A typical setup using a PhotoniQ, four photomultiplier tubes, optics, a laser, and a microcapillary electrophoresis array containing the DNA fragments is shown below. The PMTs are positioned with the optics to detect the fluorescence from the DNA fragments labeled with the individual dye markers. Each PMT connects to a charge sensitive preamp input on a PhotoniQ IQSP418 or IQSP518 multichannel PMT data acquisition system. The data acquisition system is triggered to coincide with the firing of the excitation laser such that each event consisting of the integrated outputs from the four PMTs, is digitized and recorded by the unit. The resulting data from the PhotoniQ is sent to a PC over a USB 2.0 connection for display, logging, or real time processing. Bias to the PMT detectors is controlled by connecting the front panel DAC output from the PhotoniQ to the detector bias adjust input on the PMTs. Since many PMTs have a low voltage control input to set the high voltage to the cathode, this provides a convenient means for the user to set the PMT gain through the PhotoniQ GUI or user's custom software application.



In DNA sequencing the PhotoniQ is configured in edge trigger mode so that the unit captures the fluorescence signals detected by the PMTs for a pre-defined period of time immediately following the firing of the excitation laser. This period is called the *integration time* and corresponds to the interval over which the PhotoniQ integrates the resulting fluorescence signal. Timing for this mode is shown below.





Vertilon Corporation has made every attempt to ensure that the information in this document is accurate and complete. Vertilon assumes no liability for errors or for any incidental, consequential, indirect, or special damages including, without limitation, loss of use, loss or alteration of data, delays, lost profits or savings, arising from the use of this document or the product which it accompanies.

Vertilon reserves the right to change this product without prior notice. No responsibility is assumed by Vertilon for any infringements of patents or other rights of third parties which may result from its use. No license is granted by implication or otherwise under the patent and proprietary information rights of Vertilon Corporation.

© 2010 Vertilon Corporation, ALL RIGHTS RESERVED

No form of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose without prior, express written consent from Vertilon Corporation.

AN3322.1.0 Jan 2010